


# UbiQ

targeting the ubiquitin system



## Cy5-Ub-VME (human sequence, synthetic)

UbiQ code : UbiQ-071

Batch # : B01082013-001

Amount : 50 ug, lyophilized powder

Purity :  $\geq 95\%$  by RP-HPLC

Mol. Weight : 9.1 kDa


Storage : upon arrival powder at  $-20^{\circ}\text{C}$ ; buffered solution at  $-80^{\circ}\text{C}$ . Please protect from light and avoid multiple freeze/thaw cycles.

## Productsheet

**Background.** Cy5-Ub-VME (**UbiQ-071**) is a potent, irreversible and specific inhibitor of deubiquitinating enzymes (DUBs),<sup>1</sup> which is labeled on the N-terminus with a Cy5 dye (Cy5, exc 625-650 nm, abs 670 nm).<sup>2</sup> This ubiquitin-based activity probe can be used for activity profiling experiments and the control of DUB inhibitor specificity.<sup>2,3</sup> Whereas the first-generation activity probes (such as HA tagged Ub-VME) require immunoblotting for detection, **UbiQ-071** allows for detection of DUB labeling by in-gel fluorescence. This direct and more sensitive read-out gives more distinct labeling patterns than can be achieved by immunoblotting. In addition, cross-reactivity of antibodies can lead to background labeling, something that is not observed with **UbiQ-071**.

### Sequence

Cy5-MQIFVKLTGKITLEVEPSDTIENVKAKIQDKEGIPDPQRLIFAGKQLEDGRTLSDYNIQKESTLHLVLRG-VME



**Figure 1. A; LC-MS analysis.** Mobile phase A= 1%  $\text{CH}_3\text{CN}$ , 0.1% formic acid in milliQ and B= 1% milliQ and 0.1% formic acid in  $\text{CH}_3\text{CN}$ . Phenomenex Kinetex C18, (2.1 $\times$ 50 mm, 2.6  $\mu\text{M}$ ); flow rate = 0.5 mL/min, column T =  $40^{\circ}\text{C}$ . Gradient: 5%  $\Rightarrow$  95% over 3.5 min. **B: SDS-PAGE analysis**, 12% gel, MES buffer. Left: fluorescence scanning (650/690 nm), right: CBB staining.

## Important: sample preparation

- dissolve the powder in as little DMSO as possible (e.g. 20 mg/mL)
- add this DMSO stock slowly to milliQ (please note the order of addition).
- buffer as desired (with e.g. 1M HEPES to 50 mM HEPES).
- a final buffered stock of for example 0.5 mg/mL contains 2.5 vol% DMSO; in general DMSO concentrations of up to 5 vol% are well tolerated by DUBs.
- if required, total removal of DMSO is accomplished by dialysis or spin-filtration (3 kDa cut-off membrane).

**Literature.** (1) Misaghi et al. *J. Biol. Chem.* **2005**, *280*, 1512. (2) de Jong et al. *ChemBioChem* **2012**, *13*, 2251. (3) Altun et al. *Chem. Biol.* **2011**, *18*, 1401.