

UbiQ

targeting the ubiquitin system

Biotin-Ahx-Ub-PA (human sequence, synthetic)


UbiQ code : UbiQ-076
Batch # : B01082013-001
Amount : 50 ug, lyophilized powder
Purity : $\geq 95\%$ by RP-HPLC and SDS-PAGE
Mol. Weight : 8.89 kDa
Storage : upon arrival, powder at -20°C ; buffered solution at -80°C . Please avoid multiple freeze/thaw cycles.

Productsheet

Background. UbiQ-076 is a potent and specific inhibitor of deubiquitinating enzymes (DUBs), which is labeled on the N-terminus with biotin and a C-terminal propargylamide (PA, alternative abbreviation= Prg). An aminohexanoic acid (Ahx) linker is used to create extra space between the biotin and Ub protein for efficient access of biotin binding entities. UbiQ-076 can be used for activity profiling experiments and determining DUB inhibitor specificity.¹⁻⁴ The PA part has two unique capabilities. First, it forms a covalent linkage with (the active site Cys residue of) a DUB that can be cleaved by acid treatment (5% aq. TFA), allowing for proteomic analyses. Secondly, it targets the three major DUB families: UCH, USP and OTU.^{1,2}

sequence

Biotin-Ahx-MQIFVKTLTGKTTITLEVEPSDTIENVKAKIQDKEGIPPDQQLIFAGKQLEDGRTLSDYNIQKESTLHLVLRRLRG-PA



A: LC-MS analysis. Mobile phase A= 1% CH_3CN , 0.1% formic acid in milliQ, B= 1% milliQ and 0.1% formic acid in CH_3CN . Phenomenex Kinetex C18, (2.1 \times 50 mm, 2.6 μM); flow rate= 0.5 mL/min, column T= 40°C . Gradient: 5% \Rightarrow 95% over 3.5 min. **B: SDS-PAGE analysis.** Coomassie blue staining, 12% SDS-PAGE gel.

important: sample preparation

- dissolve the powder in as little DMSO as possible (e.g. 20 mg/mL) and add this DMSO stock slowly to milliQ (please note the order of addition).
- next, buffer as desired.¹⁻³ In general, HEPES and Tris buffers are standard for DUB assays. Please note that certain DUBs react different to low or high NaCl concentrations.
- a final buffered stock of for example 0.5 mg/mL contains 2.5 vol% DMSO.
- if required, total removal of DMSO is accomplished by dialysis or spin-filtration (3 kDa cut-off membrane).
- for detailed experimental conditions please see open-access reference 1

* please note: Ub proteins tend to show aggregation during SDS-PAGE analysis of (di)Ub conjugates. This can be caused by (heat-induced) aggregation (Morimoto et al. *Sci Rep* **2018**, 8, article 2711). If possible, avoid heating the samples in Laemmli sample buffer and/or add 4M urea to your SDS-PAGE analysis samples.

Literature. (1) Ekkebus et al. *J Am Chem Soc* **2013**, 135, 2867. (2) Sommer et al. *Bioorg Med Chem* **2013**, 21, 2511. (3) de Jong et al. *ChemBioChem* **2012**, 13, 2251.