

Cy5-Ub-Dha

UbiQ code : UbiQ-104
 Batch # : B01072015-001
 Amount : 50 ug, lyophilized powder
 Purity : ≥95% by RP-HPLC
 Mol. Weight : 9.05 kDa
 Storage : powder at -20°C; solution at -80°C. Please avoid multiple freeze/thaw cycles.

Productsheet

Background. UbiQ-104 is a new and first of its kind activity based probe for Ub E1, E2 and (HECT/RBR) E3 ligases.¹ It is based on the Ub sequence in which the C-terminal Gly76 has been replaced by a dehydroalanine (Dha) residue. The N-terminus is labeled with a Cy5 tag ($\lambda_{ex} = 625 \text{ nm}$; $\lambda_{em} = 680 \text{ nm}$). It has been prepared by total chemical synthesis and is therefore well-defined in terms of Cy5 site.

UbiQ-104 is processed in a native manner by Ub E1, E2 and (HECT/RBR) E3 ligases and during this process it forms an electrophilic intermediate that can react with the active site Cys residue of the E1, E2 and (HECT/RBR) E3 enzyme, thereby creating a covalent bond (Figure 1).

Sequence

Cys-MQIFVKLTGKTITLEVPSDTIENVKAKIQDKEGIPPDQQLIFAGKQLEDGRTLSDYNIQKESTLHLVLRRLRG-Dha




Figure 1 - Mode of action of Ub-Dha activity based probes for E1-E2 and (HECT/RBR)-E3 enzymes.

Important: sample preparation

- dissolve the powder in as little DMSO as possible (e.g. 20 mg/mL)
- add this DMSO stock slowly to milliQ (please note the order of addition)
- buffer the aq. solution as desired (e.g. 50 mM HEPES pH 8, 100 mM NaCl)
- final stocks of e.g. 0.5 mg/mL will contain 2.5 vol% DMSO.
- buffer exchange using 3 kDa spin filters or dialysis membrane allows total removal of DMSO if desired; this is however not required as in general <5 vol% DMSO is well tolerated by most enzymes.

General Experimental Conditions E1 labeling assay.

UBE1 or UBA6 (1 μ M) in 50 mM HEPES pH 8, 100 mM NaCl, 10 mM MgCl₂ and 250 μ M ATP was incubated with probe (30 μ M) at 37°C for >30 min. The reaction was quenched by the addition of reducing sample buffer and heating (90°C for 10 min).

General Experimental Conditions E1+E2 labeling assay.

E2 enzyme (2.5 μ M) and UBE1 (0.63 μ M) in 50 mM HEPES pH 7.5, 100 mM NaCl, 5 mM MgCl₂ and 250 μ M ATP were incubated with probe (12.5 μ M) at 37°C for >30 min. The reaction was quenched by the addition of reducing sample buffer and heating (90°C for 10 min).

General Experimental Conditions E1+E2+(HECT/RBR) E3 labeling assay.

Nedd4L (2.5 μ M), UBE2D (0.5 μ M) and UBE1 (0.25 μ M) were incubated with probe (50 μ M) in 50 mM HEPES pH 7.5, 100 mM NaCl, 5 mM MgCl₂ and 250 μ M ATP at 30°C for 2h. The reaction was quenched by the addition of reducing sample buffer and heating (90°C for 10 min).

Cell lysate labeling.

Cell lysates are resuspended in 3 pellet volumes of HR buffer (50 mM TRIS, pH 7.4, 5 mM MgCl₂, 250 mM sucrose, 1 mM DTT) and lysed by sonication. After clarification by centrifugation (20,000 rpm, 4°C, 20 min), total protein concentration is determined by Nanodrop. For labeling experiments, 100 μ g of lysate is incubated with 0.5 μ g **Cy5-Ub-Dha**, 10 mM ATP, 10 mM MgCl₂, in Labeling Buffer (50 mM HEPES, 100 mM NaCl, pH 7.5) at 37°C for 1h or as indicated. Additional 1 mM ATP and MgCl₂ were added to the reaction every 20 minutes to replenish consumed ATP. In case of a negative control, lysates were treated with 2 Units of Apyrase (Sigma Aldrich) prior to addition of the probe (ATP and MgCl₂ were omitted). The reaction was terminated by the addition of 3x SDS-PAGE Loading Buffer (Invitrogen) containing beta-mercaptoethanol. Samples were resolved by SDS-PAGE and visualized by fluorescence scanning (λ_{ex} = 625 nm; λ_{em} = 680 nm).

Literature. (1) (a) Mulder et al. Mulder et al. *Nat. Chem. Biol.* **2016**, *12*, 523. (b) MPC Mulder, F. El Oualid and H. Ovaa. Adenylation enzyme inhibitors. Application WO/2016/032332 and NL2015/050596