

# UbiQ

targeting the ubiquitin system




Figure 1. Rh110MP

## Ac-ISG15<sup>prox</sup>-Rh110MP (mouse sequence, proximal domain, synthetic)

UbiQ code : UbiQ-127

Batch # : B01042017-001

Amount : 50 ug, lyophilized powder

Purity : ≥95%

Mol. Weight : 9.69 kDa

Storage : upon arrival, powder at -20°C, solution at -80°C. Please store dark and avoid multiple freeze/thaw cycles.

## Productsheet

**Background.** Ac-ISG15<sup>prox</sup>-Rh110MP is a quenched, fluorescent substrate for ISG15 proteases based on the proximal part of mouse ISG15. Cleavage of the amide bond between the C-terminal Gly and rhodamine110 moiety releases the highly fluorescent Rh110-morpholinecarbonyl (Figure 1, Rh110MP).<sup>1</sup>

- keep the excellent properties of the classic ubiquitin-Rh110 substrate<sup>2</sup>
- with high fluorescence intensity after proteolytic cleavage

### sequence


Ac-SEPLSILVRNERGHSNIYEVFLTQTVDTLKKKVSQREQVHEDQFWLSFEGRPMEKELLGEYGLKQPQSTVIKHLRLRGG-Rh110MP

### important: sample preparation

- dissolve the powder in DMSO: DMSO stocks may range from 1 mg/mL (103 µM) to 40 mg/mL (4.12 mM)
- add the DMSO stock to milliQ (please note the order of addition) and mix
- buffer the aq. solution as desired (using 1M HEPES or 1M Tris for example)
- a final assay stock of 100 nM will contain 0.1 vol% DMSO when prepared from a 100 µM DMSO stock

### recommended filter settings for Rh110MP

- $\lambda_{ex}$  = 492 nm,  $\lambda_{emi}$  = 525 nm, bandwidth: ± 8 nm



**LC-MS analysis.** Mobile phase A= 1% CH<sub>3</sub>CN, 0.1% formic acid in milliQ and B= 1% milliQ and 0.1% formic acid in CH<sub>3</sub>CN. XBridge BEH300 C18 3.5 µm 4.6x100mm; column T= 40°C, flow= 0.8 mL/min. Gradient: 40–75% over 6.5 min.

**Literature.** (1) Basters et al. *Nat Struct Mol Biol* **2017**, *24*, 270. (2) Terentyeva et al. *Biocon Chem* **2011**, *22*, 1932. (3) Hassiepen et al. *Analyt Biochem* **2007**, *371*, 201.