

# UbiQ

targeting the ubiquitin system

## peptide degradation assay

<b>version</b>	V1, Sept-2016
<b>applicable for</b>	UbiQ Proteasome assay reagent Mca-KKVAPYPME-Dap(Dnp)-NH <sub>2</sub>
<b>reference</b>	Jastrab, J. B., et al. An adenosine triphosphate-independent proteasome activator contributes to the virulence of <i>Mycobacterium tuberculosis</i> . <i>Proc Natl Acad Sci USA</i> . <b>14</b> , E1763-E1772 (112) doi: 10.1073/pnas.1423319112

## protocol

### peptide reconstitution

- dissolve 100 µg of UbiQ-141 in 65 µl DMSO to give a 1 mM stock solution (50×)
- this can be aliquoted and stored at -20°C.
- before an experiment, an aliquot is thawed and diluted 10× in peptide degradation assay reaction buffer (50 mM Tris pH 8, 5 mM MgCl<sub>2</sub>) to make a 5× working dilution.

### peptide degradation assay

- 1 µg freshly purified *M. tuberculosis* His6-tagged 20S core particle is mixed with 3.44 µg purified His6-tagged PafE in peptide degradation assay reaction buffer (50 mM Tris pH 8, 5 mM MgCl<sub>2</sub>) in a total volume of 240 µl.
- Mixtures are incubated at 37°C for 30 min to promote complex formation, and then incubated at room temperature for 10 minutes to cool.
- Triplicate 60 µl aliquots are transferred to a 96 well plate, and 15 µl of peptide substrate (5×) is added to each well to give a final peptide concentration of 20 µM.
- Degradation rate was assessed by monitoring fluorescence generated over time at room temperature (excitation: 340 nm; emission: 405 nm).